

BLOOD COLLECTION NEEDLES PRODUCT CODES



BLOOD COLLECTION NEEDLE, TYP I

Product code	Designation	Description
NE-BC01	COLLECTprim	Blood collection needle pen type, 23G, 1" (0.6x25mm)
NE-BC02	COLLECTprim	Blood collection needle pen type, 23G, 1 1/2" (0.6x38mm)
NE-BC03	COLLECTprim	Blood collection needle pen type, 22G, 1" (0.7x25mm)
NE-BC04	COLLECTprim	Blood collection needle pen type, 22G, 1 1/2" (0.7x38mm)
NE-BC05	COLLECTprim	Blood collection needle pen type, 21G, 1" (0.8x25mm)
NE-BC06	COLLECTprim	Blood collection needle pen type, 21G, 1 1/2" (0.8x38mm)
NE-BC07	COLLECTprim	Blood collection needle pen type, 20G, 1" (0.9x25mm)
NE-BC08	COLLECTprim	Blood collection needle pen type, 20G, 1 1/2" (0.9x38mm)
NE-BC09	COLLECTprim	Blood collection needle pen type, 18G, 1" (1.2x25mm)
NE-BC10	COLLECTprim	Blood collection needle pen type, 18G, 1 1/2" (1.2x38mm)

BLOOD COLLECTION NEEDLE, TYPE II

Product code	Designation	Description
NE-BC11	COLLECTprim	Blood collection needle with scalp vein set, 24G, 3/4" (0,55x20mm), long bevel
NE-BC12	COLLECTprim	Blood collection needle with scalp vein set, 24G, 3/4" (0,55x20mm), short bevel
NE-BC13	COLLECTprim	Blood collection needle with scalp vein set, 24G, 1 1/2" (0,55x38mm), long bevel
NE-BC14	COLLECTprim	Blood collection needle with scalp vein set, 24G, 1 1/2" (0,55x38mm), short bevel
NE-BC15	COLLECTprim	Blood collection needle with scalp vein set, 23G, 3/4" (0,60x20mm), long bevel
NE-BC16	COLLECTprim	Blood collection needle with scalp vein set, 23G, 3/4" (0,60x20mm), short bevel
NE-BC17	COLLECTprim	Blood collection needle with scalp vein set, 23G, 1 1/2" (0,60x38mm), long bevel
NE-BC18	COLLECTprim	Blood collection needle with scalp vein set, 23G, 1 1/2" (0,60x38mm), short bevel
NE-BC19	COLLECTprim	Blood collection needle with scalp vein set, 22G, 3/4" (0,70x20mm), long bevel
NE-BC20	COLLECTprim	Blood collection needle with scalp vein set, 22G, 3/4" (0,70x20mm), short bevel
NE-BC21	COLLECTprim	Blood collection needle with scalp vein set, 22G, 1 1/2" (0,70x38mm), long bevel
NE-BC22	COLLECTprim	Blood collection needle with scalp vein set, 22G, 1 1/2" (0,70x38mm), short bevel
NE-BC23	COLLECTprim	Blood collection needle with scalp vein set, 21G, 3/4" (0,80x20mm), long bevel
NE-BC24	COLLECTprim	Blood collection needle with scalp vein set, 21G, 3/4" (0,80x20mm), short bevel
NE-BC25	COLLECTprim	Blood collection needle with scalp vein set, 21G, 1 1/2" (0,80x38mm), long bevel
NE-BC26	COLLECTprim	Blood collection needle with scalp vein set, 21G, 1 1/2" (0,80x38mm), short bevel
NE-BC27	COLLECTprim	Blood collection needle with scalp vein set, 20G, 3/4" (0,90x20mm), long bevel
NE-BC28	COLLECTprim	Blood collection needle with scalp vein set, 20G, 3/4" (0,90x20mm), short bevel
NE-BC29	COLLECTprim	Blood collection needle with scalp vein set, 20G, 1 1/2" (0,90x38mm), long bevel
NE-BC30	COLLECTprim	Blood collection needle with scalp vein set, 20G, 1 1/2" (0,90x38mm), short bevel
NE-BC31	COLLECTprim	Blood collection needle with scalp vein set, 19G, 3/4" (1,1x20mm), long bevel
NE-BC32	COLLECTprim	Blood collection needle with scalp vein set, 19G, 3/4" (1,1x20mm), short bevel
NE-BC33	COLLECTprim	Blood collection needle with scalp vein set, 19G, 1 1/2" (1,1x38mm), long bevel
NE-BC34	COLLECTprim	Blood collection needle with scalp vein set, 19G, 1 1/2" (1,1x38mm), short bevel
NE-BC35	COLLECTprim	Blood collection needle with scalp vein set, 18G, 3/4" (1,2x20mm), long bevel
NE-BC36	COLLECTprim	Blood collection needle with scalp vein set, 18G, 3/4" (1,2x20mm), short bevel
NE-BC37	COLLECTprim	Blood collection needle with scalp vein set, 18G, 1 1/2" (1,2x38mm), long bevel
NE-BC38	COLLECTprim	Blood collection needle with scalp vein set, 18G, 1 1/2" (1,2x38mm), short bevel