


LIAprim INFUSION PUMPS PRODUCT CODES



TYPE I

This type of LIAprim infusion pumps has a continuous flow rate (also known as Continuous Basal Infusion or CBI type).



Description

- 1 Carrying device (bracket and a fabric cord)
- 2 Plastic container with balloon reservoir
- 3 Luer-lock inlet with a cap
- 4 Tubing
- 5 Infusion filter
- 6 Luer-lock connector
- 7 Protector cap

Configurations


Volume	Flow rate				
	1ml/h	2ml/h	4ml/h	5ml/h	10ml/h
60 ml	x	x			
100 ml	x	x	x	x	x
150 ml	x	x	x	x	
200 ml	x	x	x	x	
275 ml	x	x	x	x	x
300 ml	x	x	x	x	

PACKAGING

Sterilisation: Ethylene oxide
 Shelf life: 5 years after manufacturing

TYPE II

This type of LIAprim infusion pumps has an additional flow rate regulating device providing four different flow rates to choose from.



Description

- 1 Carrying device (bracket and a fabric cord)
- 2 Plastic container with balloon reservoir
- 3 Luer-lock inlet with a cap
- 4 Tubing
- 5 Infusion filter
- 6 Flow rate regulating device
- 7 Luer-lock connector
- 8 Protector cap

Configurations


	Flow rate				
	1/2/3/ 4ml/h	2/3/4/ 5ml/h	2/4/6/ 8ml/h	4/6/8/ 10ml/h	6/8/10/ 12ml/h
Volume					
60 ml	x		x		
100 ml	x	x	x	x	
150 ml	x		x	x	
200 ml		x	x	x	x
275 ml			x	x	x
300 ml		x	x	x	x

PACKAGING

Sterilisation: Ethylene oxide
 Shelf life: 5 years after manufacturing

TYPE III

This type of LIAprim infusion pumps has an additional bolus function, so that the patient can self-administer an extra dosage of analgesic according to actual pain situation.



Description

- 1 Carrying device (bracket and a fabric cord)
- 2 Plastic container with balloon reservoir
- 3 Luer-lock inlet with a cap
- 4 Tubing
- 5 Infusion filter
- 6 Bolus
- 7 Luer-lock connector
- 8 Protector cap

Configurations

Volume	Flow rate				
	1ml/h	2ml/h	4ml/h	5ml/h	10ml/h
60 ml	x	x			
100 ml	x	x	x	x	x
150 ml	x	x	x	x	
200 ml	x	x	x	x	
275 ml	x	x	x	x	x
300 ml	x	x	x	x	

EVERY CONFIGURATION IS AVAILABLE WITH ALL 6 POSSIBLE TYPES OF BOLUS (DOSAGE / LOCK OUT TIME)


- 0,5 ml/10 min
- 0,5 ml/15 min
- 1,0 ml/10 min
- 1,0 ml/15 min
- 2,0 ml/10 min
- 2,0 ml/15 min

PACKAGING

Sterilisation: Ethylene oxide
 Shelf life: 5 years after manufacturing

TYPE IV

This type of LIAprim infusion pumps has both: a flow rate regulating device providing four different flow rates to choose from and a bolus function, so that the patient can self-administer an extra dosage of analgesic according to actual pain situation.



Description

- 1 Carrying device (bracket and a fabric cord)
- 2 Plastic container with balloon reservoir
- 3 Luer-lock inlet with a cap
- 4 Tubing
- 5 Infusion filter
- 6 Flow rate regulating device
- 7 Bolus
- 8 Luer-lock connector
- 9 Protector cap

Configurations

	Flow rate				
	1/2/3/ 4ml/h	2/3/4/ 5ml/h	2/4/6/ 8ml/h	4/6/8/ 10ml/h	6/8/10/ 12ml/h
Volume					
60 ml	x		x	x	
100 ml	x		x	x	
150 ml	x		x	x	
200 ml			x	x	
275 ml			x	x	x
300 ml			x	x	x

EVERY CONFIGURATION IS AVAILABLE WITH ALL 6 POSSIBLE TYPES OF BOLUS (DOSAGE / LOCK OUT TIME)

- 0,5 ml/10 min
- 0,5 ml/15 min
- 1,0 ml/10 min
- 1,0 ml/15 min
- 2,0 ml/10 min
- 2,0 ml/15 min

The most popular configurations of LIAprim infusion pumps, type IV you can find in the Product code table (P. 98). All other configurations available on request.

PACKAGING

- Sterilisation: Ethylene oxide
- Shelf life: 5 years after manufacturing

LIAprim INFUSION PUMP, TYPE I

Product code	Designation	Description
IS-IP01	LIAprim	Infusion pump 60 ml, flow rate 1 ml/h
IS-IP02	LIAprim	Infusion pump 60 ml, flow rate 2 ml/h
IS-IP10	LIAprim	Infusion pump 100 ml, flow rate 1 ml/h
IS-IP11	LIAprim	Infusion pump 100 ml, flow rate 2 ml/h
IS-IP13	LIAprim	Infusion pump 100 ml, flow rate 4 ml/h
IS-IP14	LIAprim	Infusion pump 100 ml, flow rate 5 ml/h
IS-IP17	LIAprim	Infusion pump 100 ml, flow rate 10 ml/h
IS-IP19	LIAprim	Infusion pump 150 ml, flow rate 1 ml/h
IS-IP20	LIAprim	Infusion pump 150 ml, flow rate 2 ml/h
IS-IP22	LIAprim	Infusion pump 150 ml, flow rate 4 ml/h
IS-IP23	LIAprim	Infusion pump 150 ml, flow rate 5 ml/h
IS-IP28	LIAprim	Infusion pump 200 ml, flow rate 1 ml/h
IS-IP29	LIAprim	Infusion pump 200 ml, flow rate 2 ml/h
IS-IP31	LIAprim	Infusion pump 200 ml, flow rate 4 ml/h
IS-IP32	LIAprim	Infusion pump 200 ml, flow rate 5 ml/h
IS-IP37	LIAprim	Infusion pump 275 ml, flow rate 1 ml/h
IS-IP38	LIAprim	Infusion pump 275 ml, flow rate 2 ml/h
IS-IP40	LIAprim	Infusion pump 275 ml, flow rate 4 ml/h
IS-IP41	LIAprim	Infusion pump 275 ml, flow rate 5 ml/h
IS-IP44	LIAprim	Infusion pump 275 ml, flow rate 10 ml/h
IS-IP46	LIAprim	Infusion pump 300 ml, flow rate 1 ml/h
IS-IP47	LIAprim	Infusion pump 300 ml, flow rate 2 ml/h
IS-IP49	LIAprim	Infusion pump 300 ml, flow rate 4 ml/h
IS-IP50	LIAprim	Infusion pump 300 ml, flow rate 5 ml/h

LIAprim INFUSION PUMP, TYPE II

Product code	Designation	Description
IS-IP381	LIAprim	Infusion pump 60 ml, dosing volume 1/2/3/4ml/h
IS-IP383	LIAprim	Infusion pump 60 ml, dosing volume 2/4/6/8ml/h
IS-IP386	LIAprim	Infusion pump 100 ml, dosing volume 1/2/3/4ml/h
IS-IP387	LIAprim	Infusion pump 100 ml, dosing volume 2/3/4/5ml/h
IS-IP388	LIAprim	Infusion pump 100 ml, dosing volume 2/4/6/8ml/h
IS-IP389	LIAprim	Infusion pump 100 ml, dosing volume 4/6/8/10ml/h
IS-IP391	LIAprim	Infusion pump 150 ml, dosing volume 1/2/3/4ml/h
IS-IP393	LIAprim	Infusion pump 150 ml, dosing volume 2/4/6/8ml/h
IS-IP394	LIAprim	Infusion pump 150 ml, dosing volume 4/6/8/10ml/h
IS-IP397	LIAprim	Infusion pump 200 ml, dosing volume 2/3/4/5ml/h
IS-IP398	LIAprim	Infusion pump 200 ml, dosing volume 2/4/6/8ml/h
IS-IP399	LIAprim	Infusion pump 200 ml, dosing volume 4/6/8/10ml/h
IS-IP400	LIAprim	Infusion pump 200 ml, dosing volume 6/8/10/12ml/h
IS-IP403	LIAprim	Infusion pump 275 ml, dosing volume 2/4/6/8ml/h
IS-IP404	LIAprim	Infusion pump 275 ml, dosing volume 4/6/8/10ml/h
IS-IP404	LIAprim	Infusion pump 275 ml, dosing volume 6/8/10/12ml/h
IS-IP407	LIAprim	Infusion pump 300 ml, dosing volume 2/3/4/5ml/h
IS-IP408	LIAprim	Infusion pump 300 ml, dosing volume 2/4/6/8ml/h
IS-IP409	LIAprim	Infusion pump 300 ml, dosing volume 4/6/8/10ml/h
IS-IP410	LIAprim	Infusion pump 300 ml, dosing volume 6/8/10/12ml/h

LIAprim INFUSION PUMP, TYPE III

Product code	Designation	Description
IS-IP55	LIAprim	Infusion pump 60 ml, flow rate 1 ml/h, PCA Bolus 0.5ml/10min
IS-IP56	LIAprim	Infusion pump 60 ml, flow rate 1 ml/h, PCA Bolus 0.5ml/15min
IS-IP57	LIAprim	Infusion pump 60 ml, flow rate 2 ml/h, PCA Bolus 0.5ml/10min
IS-IP58	LIAprim	Infusion pump 60 ml, flow rate 2 ml/h, PCA Bolus 0.5ml/15min
IS-IP165	LIAprim	Infusion pump 60 ml, flow rate 1 ml/h, PCA Bolus 1ml/10min
IS-IP166	LIAprim	Infusion pump 60 ml, flow rate 1 ml/h, PCA Bolus 1ml/15min
IS-IP167	LIAprim	Infusion pump 60 ml, flow rate 2 ml/h, PCA Bolus 1ml/10min
IS-IP168	LIAprim	Infusion pump 60 ml, flow rate 2 ml/h, PCA Bolus 1ml/15min
IS-IP273	LIAprim	Infusion pump 60 ml, flow rate 1 ml/h, PCA Bolus 2ml/10min
IS-IP274	LIAprim	Infusion pump 60 ml, flow rate 1 ml/h, PCA Bolus 2ml/15min
IS-IP275	LIAprim	Infusion pump 60 ml, flow rate 2 ml/h, PCA Bolus 2ml/10min
IS-IP276	LIAprim	Infusion pump 60 ml, flow rate 2 ml/h, PCA Bolus 2ml/15min
IS-IP74	LIAprim	Infusion pump 100 ml, flow rate 1 ml/h, PCA Bolus 0.5ml/10min
IS-IP75	LIAprim	Infusion pump 100 ml, flow rate 1 ml/h, PCA Bolus 0.5ml/15min
IS-IP76	LIAprim	Infusion pump 100 ml, flow rate 2 ml/h, PCA Bolus 0.5ml/10min
IS-IP77	LIAprim	Infusion pump 100 ml, flow rate 2 ml/h, PCA Bolus 0.5ml/15min
IS-IP80	LIAprim	Infusion pump 100 ml, flow rate 4 ml/h, PCA Bolus 0.5ml/10min
IS-IP81	LIAprim	Infusion pump 100 ml, flow rate 4 ml/h, PCA Bolus 0.5ml/15min
IS-IP82	LIAprim	Infusion pump 100 ml, flow rate 5 ml/h, PCA Bolus 0.5ml/10min
IS-IP83	LIAprim	Infusion pump 100 ml, flow rate 5 ml/h, PCA Bolus 0.5ml/15min
IS-IP88	LIAprim	Infusion pump 100 ml, flow rate 10 ml/h, PCA Bolus 0.5ml/10min
IS-IP89	LIAprim	Infusion pump 100 ml, flow rate 10 ml/h, PCA Bolus 0.5ml/15min
IS-IP183	LIAprim	Infusion pump 100 ml, flow rate 1 ml/h, PCA Bolus 1ml/10min
IS-IP184	LIAprim	Infusion pump 100 ml, flow rate 1 ml/h, PCA Bolus 1ml/15min
IS-IP185	LIAprim	Infusion pump 100 ml, flow rate 2 ml/h, PCA Bolus 1ml/10min
IS-IP186	LIAprim	Infusion pump 100 ml, flow rate 2 ml/h, PCA Bolus 1ml/15min
IS-IP189	LIAprim	Infusion pump 100 ml, flow rate 4 ml/h, PCA Bolus 1ml/10min
IS-IP190	LIAprim	Infusion pump 100 ml, flow rate 4 ml/h, PCA Bolus 1ml/15min
IS-IP191	LIAprim	Infusion pump 100 ml, flow rate 5 ml/h, PCA Bolus 1ml/10min
IS-IP192	LIAprim	Infusion pump 100 ml, flow rate 5 ml/h, PCA Bolus 1ml/15min

LIAprim INFUSION PUMP, TYPE III

Product code	Designation	Description
IS-IP197	LIAprim	Infusion pump 100 ml, flow rate 10 ml/h, PCA Bolus 1ml/10min
IS-IP198	LIAprim	Infusion pump 100 ml, flow rate 10 ml/h, PCA Bolus 1ml/15min
IS-IP291	LIAprim	Infusion pump 100 ml, flow rate 1 ml/h, PCA Bolus 2ml/10min
IS-IP292	LIAprim	Infusion pump 100 ml, flow rate 1 ml/h, PCA Bolus 2ml/15min
IS-IP293	LIAprim	Infusion pump 100 ml, flow rate 2 ml/h, PCA Bolus 2ml/10min
IS-IP294	LIAprim	Infusion pump 100 ml, flow rate 2 ml/h, PCA Bolus 2ml/15min
IS-IP297	LIAprim	Infusion pump 100 ml, flow rate 4 ml/h, PCA Bolus 2ml/10min
IS-IP298	LIAprim	Infusion pump 100 ml, flow rate 4 ml/h, PCA Bolus 2ml/15min
IS-IP299	LIAprim	Infusion pump 100 ml, flow rate 5 ml/h, PCA Bolus 2ml/10min
IS-IP300	LIAprim	Infusion pump 100 ml, flow rate 5 ml/h, PCA Bolus 2ml/15min
IS-IP305	LIAprim	Infusion pump 100 ml, flow rate 10 ml/h, PCA Bolus 2ml/10min
IS-IP306	LIAprim	Infusion pump 100 ml, flow rate 10 ml/h, PCA Bolus 2ml/15min
IS-IP92	LIAprim	Infusion pump 150 ml, flow rate 1 ml/h, PCA Bolus 0.5ml/10min
IS-IP93	LIAprim	Infusion pump 150 ml, flow rate 1 ml/h, PCA Bolus 0.5ml/15min
IS-IP94	LIAprim	Infusion pump 150 ml, flow rate 2 ml/h, PCA Bolus 0.5ml/10min
IS-IP95	LIAprim	Infusion pump 150 ml, flow rate 2 ml/h, PCA Bolus 0.5ml/15min
IS-IP98	LIAprim	Infusion pump 150 ml, flow rate 4 ml/h, PCA Bolus 0.5ml/10min
IS-IP99	LIAprim	Infusion pump 150 ml, flow rate 4 ml/h, PCA Bolus 0.5ml/15min
IS-IP100	LIAprim	Infusion pump 150 ml, flow rate 5 ml/h, PCA Bolus 0.5ml/10min
IS-IP101	LIAprim	Infusion pump 150 ml, flow rate 5 ml/h, PCA Bolus 0.5ml/15min
IS-IP201	LIAprim	Infusion pump 150 ml, flow rate 1 ml/h, PCA Bolus 1ml/10min
IS-IP202	LIAprim	Infusion pump 150 ml, flow rate 1 ml/h, PCA Bolus 1ml/15min
IS-IP203	LIAprim	Infusion pump 150 ml, flow rate 2 ml/h, PCA Bolus 1ml/10min
IS-IP204	LIAprim	Infusion pump 150 ml, flow rate 2 ml/h, PCA Bolus 1ml/15min
IS-IP207	LIAprim	Infusion pump 150 ml, flow rate 4 ml/h, PCA Bolus 1ml/10min
IS-IP208	LIAprim	Infusion pump 150 ml, flow rate 4 ml/h, PCA Bolus 1ml/15min
IS-IP209	LIAprim	Infusion pump 150 ml, flow rate 5 ml/h, PCA Bolus 1ml/10min
IS-IP210	LIAprim	Infusion pump 150 ml, flow rate 5 ml/h, PCA Bolus 1ml/15min
IS-IP309	LIAprim	Infusion pump 150 ml, flow rate 1 ml/h, PCA Bolus 2ml/10min
IS-IP310	LIAprim	Infusion pump 150 ml, flow rate 1 ml/h, PCA Bolus 2ml/15min

LIAprim INFUSION PUMP, TYPE III

Product code	Designation	Description
IS-IP311	LIAprim	Infusion pump 150 ml, flow rate 2 ml/h, PCA Bolus 2ml/10min
IS-IP312	LIAprim	Infusion pump 150 ml, flow rate 2 ml/h, PCA Bolus 2ml/15min
IS-IP315	LIAprim	Infusion pump 150 ml, flow rate 4 ml/h, PCA Bolus 2ml/10min
IS-IP316	LIAprim	Infusion pump 150 ml, flow rate 4 ml/h, PCA Bolus 2ml/15min
IS-IP317	LIAprim	Infusion pump 150 ml, flow rate 5 ml/h, PCA Bolus 2ml/10min
IS-IP318	LIAprim	Infusion pump 150 ml, flow rate 5 ml/h, PCA Bolus 2ml/15min
IS-IP110	LIAprim	Infusion pump 200 ml, flow rate 1 ml/h, PCA Bolus 0.5ml/10min
IS-IP111	LIAprim	Infusion pump 200 ml, flow rate 1 ml/h, PCA Bolus 0.5ml/15min
IS-IP112	LIAprim	Infusion pump 200 ml, flow rate 2 ml/h, PCA Bolus 0.5ml/10min
IS-IP113	LIAprim	Infusion pump 200 ml, flow rate 2 ml/h, PCA Bolus 0.5ml/15min
IS-IP116	LIAprim	Infusion pump 200 ml, flow rate 4 ml/h, PCA Bolus 0.5ml/10min
IS-IP117	LIAprim	Infusion pump 200 ml, flow rate 4 ml/h, PCA Bolus 0.5ml/15min
IS-IP118	LIAprim	Infusion pump 200 ml, flow rate 5 ml/h, PCA Bolus 0.5ml/10min
IS-IP119	LIAprim	Infusion pump 200 ml, flow rate 5 ml/h, PCA Bolus 0.5ml/15min
IS-IP219	LIAprim	Infusion pump 200 ml, flow rate 1 ml/h, PCA Bolus 1ml/10min
IS-IP220	LIAprim	Infusion pump 200 ml, flow rate 1 ml/h, PCA Bolus 1ml/15min
IS-IP221	LIAprim	Infusion pump 200 ml, flow rate 2 ml/h, PCA Bolus 1ml/10min
IS-IP222	LIAprim	Infusion pump 200 ml, flow rate 2 ml/h, PCA Bolus 1ml/15min
IS-IP225	LIAprim	Infusion pump 200 ml, flow rate 4 ml/h, PCA Bolus 1ml/10min
IS-IP226	LIAprim	Infusion pump 200 ml, flow rate 4 ml/h, PCA Bolus 1ml/15min
IS-IP227	LIAprim	Infusion pump 200 ml, flow rate 5 ml/h, PCA Bolus 1ml/10min
IS-IP228	LIAprim	Infusion pump 200 ml, flow rate 5 ml/h, PCA Bolus 1ml/15min
IS-IP327	LIAprim	Infusion pump 200 ml, flow rate 1 ml/h, PCA Bolus 2ml/10min
IS-IP328	LIAprim	Infusion pump 200 ml, flow rate 1 ml/h, PCA Bolus 2ml/15min
IS-IP329	LIAprim	Infusion pump 200 ml, flow rate 2 ml/h, PCA Bolus 2ml/10min
IS-IP330	LIAprim	Infusion pump 200 ml, flow rate 2 ml/h, PCA Bolus 2ml/15min
IS-IP333	LIAprim	Infusion pump 200 ml, flow rate 4 ml/h, PCA Bolus 2ml/10min
IS-IP334	LIAprim	Infusion pump 200 ml, flow rate 4 ml/h, PCA Bolus 2ml/15min
IS-IP335	LIAprim	Infusion pump 200 ml, flow rate 5 ml/h, PCA Bolus 2ml/10min
IS-IP336	LIAprim	Infusion pump 200 ml, flow rate 5 ml/h, PCA Bolus 2ml/15min

LIAprim INFUSION PUMP, TYPE III

Product code	Designation	Description
IS-IP128	LIAprim	Infusion pump 275 ml, flow rate 1 ml/h, PCA Bolus 0.5ml/10min
IS-IP129	LIAprim	Infusion pump 275 ml, flow rate 1 ml/h, PCA Bolus 0.5ml/15min
IS-IP130	LIAprim	Infusion pump 275 ml, flow rate 2 ml/h, PCA Bolus 0.5ml/10min
IS-IP131	LIAprim	Infusion pump 275 ml, flow rate 2 ml/h, PCA Bolus 0.5ml/15min
IS-IP134	LIAprim	Infusion pump 275 ml, flow rate 4 ml/h, PCA Bolus 0.5ml/10min
IS-IP135	LIAprim	Infusion pump 275 ml, flow rate 4 ml/h, PCA Bolus 0.5ml/15min
IS-IP136	LIAprim	Infusion pump 275 ml, flow rate 5 ml/h, PCA Bolus 0.5ml/10min
IS-IP137	LIAprim	Infusion pump 275 ml, flow rate 5 ml/h, PCA Bolus 0.5ml/15min
IS-IP142	LIAprim	Infusion pump 275 ml, flow rate 10 ml/h, PCA Bolus 0.5ml/10min
IS-IP143	LIAprim	Infusion pump 275 ml, flow rate 10 ml/h, PCA Bolus 0.5ml/15min
IS-IP237	LIAprim	Infusion pump 275 ml, flow rate 1 ml/h, PCA Bolus 1ml/10min
IS-IP238	LIAprim	Infusion pump 275 ml, flow rate 1 ml/h, PCA Bolus 1ml/15min
IS-IP239	LIAprim	Infusion pump 275 ml, flow rate 2 ml/h, PCA Bolus 1ml/10min
IS-IP240	LIAprim	Infusion pump 275 ml, flow rate 2 ml/h, PCA Bolus 1ml/15min
IS-IP243	LIAprim	Infusion pump 275 ml, flow rate 4 ml/h, PCA Bolus 1ml/10min
IS-IP244	LIAprim	Infusion pump 275 ml, flow rate 4 ml/h, PCA Bolus 1ml/15min
IS-IP245	LIAprim	Infusion pump 275 ml, flow rate 5 ml/h, PCA Bolus 1ml/10min
IS-IP246	LIAprim	Infusion pump 275 ml, flow rate 5 ml/h, PCA Bolus 1ml/15min
IS-IP251	LIAprim	Infusion pump 275 ml, flow rate 10 ml/h, PCA Bolus 1ml/10min
IS-IP252	LIAprim	Infusion pump 275 ml, flow rate 10 ml/h, PCA Bolus 1ml/15min
IS-IP345	LIAprim	Infusion pump 275 ml, flow rate 1 ml/h, PCA Bolus 2ml/10min
IS-IP346	LIAprim	Infusion pump 275 ml, flow rate 1 ml/h, PCA Bolus 2ml/15min
IS-IP347	LIAprim	Infusion pump 275 ml, flow rate 2 ml/h, PCA Bolus 2ml/10min
IS-IP348	LIAprim	Infusion pump 275 ml, flow rate 2 ml/h, PCA Bolus 2ml/15min
IS-IP351	LIAprim	Infusion pump 275 ml, flow rate 4 ml/h, PCA Bolus 2ml/10min
IS-IP352	LIAprim	Infusion pump 275 ml, flow rate 4 ml/h, PCA Bolus 2ml/15min
IS-IP353	LIAprim	Infusion pump 275 ml, flow rate 5 ml/h, PCA Bolus 2ml/10min
IS-IP354	LIAprim	Infusion pump 275 ml, flow rate 5 ml/h, PCA Bolus 2ml/15min
IS-IP359	LIAprim	Infusion pump 275 ml, flow rate 10 ml/h, PCA Bolus 2ml/10min
IS-IP360	LIAprim	Infusion pump 275 ml, flow rate 10 ml/h, PCA Bolus 2ml/15min

LIAprim INFUSION PUMP, TYPE III

Product code	Designation	Description
IS-IP146	LIAprim	Infusion pump 300 ml, flow rate 1 ml/h, PCA Bolus 0.5ml/10min
IS-IP147	LIAprim	Infusion pump 300 ml, flow rate 1 ml/h, PCA Bolus 0.5ml/15min
IS-IP148	LIAprim	Infusion pump 300 ml, flow rate 2 ml/h, PCA Bolus 0.5ml/10min
IS-IP149	LIAprim	Infusion pump 300 ml, flow rate 2 ml/h, PCA Bolus 0.5ml/15min
IS-IP153	LIAprim	Infusion pump 300 ml, flow rate 4 ml/h, PCA Bolus 0.5ml/10min
IS-IP154	LIAprim	Infusion pump 300 ml, flow rate 4 ml/h, PCA Bolus 0.5ml/15min
IS-IP155	LIAprim	Infusion pump 300 ml, flow rate 5 ml/h, PCA Bolus 0.5ml/10min
IS-IP156	LIAprim	Infusion pump 300 ml, flow rate 5 ml/h, PCA Bolus 0.5ml/15min
IS-IP255	LIAprim	Infusion pump 300 ml, flow rate 1 ml/h, PCA Bolus 1ml/10min
IS-IP256	LIAprim	Infusion pump 300 ml, flow rate 1 ml/h, PCA Bolus 1ml/15min
IS-IP257	LIAprim	Infusion pump 300 ml, flow rate 2 ml/h, PCA Bolus 1ml/10min
IS-IP258	LIAprim	Infusion pump 300 ml, flow rate 2 ml/h, PCA Bolus 1ml/15min
IS-IP261	LIAprim	Infusion pump 300 ml, flow rate 4 ml/h, PCA Bolus 1ml/10min
IS-IP262	LIAprim	Infusion pump 300 ml, flow rate 4 ml/h, PCA Bolus 1ml/15min
IS-IP263	LIAprim	Infusion pump 300 ml, flow rate 5 ml/h, PCA Bolus 1ml/10min
IS-IP264	LIAprim	Infusion pump 300 ml, flow rate 5 ml/h, PCA Bolus 1ml/15min
IS-IP363	LIAprim	Infusion pump 300 ml, flow rate 1 ml/h, PCA Bolus 2ml/10min
IS-IP364	LIAprim	Infusion pump 300 ml, flow rate 1 ml/h, PCA Bolus 2ml/15min
IS-IP365	LIAprim	Infusion pump 300 ml, flow rate 2 ml/h, PCA Bolus 2ml/10min
IS-IP366	LIAprim	Infusion pump 300 ml, flow rate 2 ml/h, PCA Bolus 2ml/15min
IS-IP369	LIAprim	Infusion pump 300 ml, flow rate 4 ml/h, PCA Bolus 2ml/10min
IS-IP370	LIAprim	Infusion pump 300 ml, flow rate 4 ml/h, PCA Bolus 2ml/15min
IS-IP371	LIAprim	Infusion pump 300 ml, flow rate 5 ml/h, PCA Bolus 2ml/10min
IS-IP372	LIAprim	Infusion pump 300 ml, flow rate 5 ml/h, PCA Bolus 2ml/15min

LIAprim INFUSION PUMP, TYPE IV

Product code	Designation	Description
IS-IP680	LIAprim	Infusion pump 60ml, flow rate 1/2/3/4ml/h, PCA Bolus 0,5ml/10min
IS-IP682	LIAprim	Infusion pump 60ml, flow rate 1/2/3/4ml/h, PCA Bolus 1ml/10min
IS-IP692	LIAprim	Infusion pump 60ml, flow rate 2/4/6/8ml/h, PCA Bolus 0,5ml/10min
IS-IP694	LIAprim	Infusion pump 60ml, flow rate 2/4/6/8ml/h, PCA Bolus 1ml/10min
IS-IP698	LIAprim	Infusion pump 60ml, flow rate 4/6/8/10ml/h, PCA Bolus 0,5ml/10min
IS-IP700	LIAprim	Infusion pump 60ml, flow rate 4/6/8/10ml/h, PCA Bolus 1ml/10min
IS-IP710	LIAprim	Infusion pump 100ml, flow rate 1/2/3/4ml/h, PCA Bolus 0,5ml/10min
IS-IP712	LIAprim	Infusion pump 100ml, flow rate 1/2/3/4ml/h, PCA Bolus 1ml/10min
IS-IP714	LIAprim	Infusion pump 100ml, flow rate 1/2/3/4ml/h, PCA Bolus 2ml/10min
IS-IP722	LIAprim	Infusion pump 100ml, flow rate 2/4/6/8ml/h, PCA Bolus 0,5ml/10min
IS-IP724	LIAprim	Infusion pump 100ml, flow rate 2/4/6/8ml/h, PCA Bolus 1ml/10min
IS-IP726	LIAprim	Infusion pump 100ml, flow rate 2/4/6/8ml/h, PCA Bolus 2ml/10min
IS-IP728	LIAprim	Infusion pump 100ml, flow rate 4/6/8/10ml/h, PCA Bolus 0,5ml/10min
IS-IP730	LIAprim	Infusion pump 100ml, flow rate 4/6/8/10ml/h, PCA Bolus 1ml/10min
IS-IP732	LIAprim	Infusion pump 100ml, flow rate 4/6/8/10ml/h, PCA Bolus 2ml/10min
IS-IP742	LIAprim	Infusion pump 150ml, flow rate 1/2/3/4ml/h, PCA Bolus 1ml/10min
IS-IP744	LIAprim	Infusion pump 150ml, flow rate 1/2/3/4ml/h, PCA Bolus 2ml/10min
IS-IP754	LIAprim	Infusion pump 150ml, flow rate 2/4/6/8ml/h, PCA Bolus 1ml/10min
IS-IP756	LIAprim	Infusion pump 150ml, flow rate 2/4/6/8ml/h, PCA Bolus 2ml/10min
IS-IP760	LIAprim	Infusion pump 150ml, flow rate 4/6/8/10ml/h, PCA Bolus 1ml/10min
IS-IP762	LIAprim	Infusion pump 150ml, flow rate 4/6/8/10ml/h, PCA Bolus 2ml/10min
IS-IP784	LIAprim	Infusion pump 200ml, flow rate 2/4/6/8ml/h, PCA Bolus 1ml/10min
IS-IP786	LIAprim	Infusion pump 200ml, flow rate 2/4/6/8ml/h, PCA Bolus 2ml/10min
IS-IP790	LIAprim	Infusion pump 200ml, flow rate 4/6/8/10ml/h, PCA Bolus 1ml/10min
IS-IP792	LIAprim	Infusion pump 200ml, flow rate 4/6/8/10ml/h, PCA Bolus 2ml/10min
IS-IP796	LIAprim	Infusion pump 200ml, flow rate 6/8/10/12ml/h, PCA Bolus 1ml/10min
IS-IP798	LIAprim	Infusion pump 200ml, flow rate 6/8/10/12ml/h, PCA Bolus 2ml/10min
IS-IP814	LIAprim	Infusion pump 275ml, flow rate 2/4/6/8ml/h, PCA Bolus 1ml/10min
IS-IP152	LIAprim	Infusion pump 275ml, flow rate 2/4/6/8ml/h, PCA Bolus 1ml/15min
IS-IP815	LIAprim	Infusion pump 275ml, flow rate 2/4/6/8ml/h, PCA Bolus 2ml/10min

LIAprim INFUSION PUMP, TYPE IV

Product code	Designation	Description
IS-IP816	LIAprim	Infusion pump 275ml, flow rate 2/4/6/8ml/h, PCA Bolus 2ml/15min
IS-IP819	LIAprim	Infusion pump 275ml, flow rate 4/6/8/10ml/h, PCA Bolus 1ml/10min
IS-IP820	LIAprim	Infusion pump 275ml, flow rate 4/6/8/10ml/h, PCA Bolus 1ml/15min
IS-IP821	LIAprim	Infusion pump 275ml, flow rate 4/6/8/10ml/h, PCA Bolus 2ml/10min
IS-IP822	LIAprim	Infusion pump 275ml, flow rate 4/6/8/10ml/h, PCA Bolus 2ml/15min
IS-IP825	LIAprim	Infusion pump 275ml, flow rate 6/8/10/12ml/h, PCA Bolus 1ml/10min
IS-IP826	LIAprim	Infusion pump 275ml, flow rate 6/8/10/12ml/h, PCA Bolus 1ml/15min
IS-IP827	LIAprim	Infusion pump 275ml, flow rate 6/8/10/12ml/h, PCA Bolus 2ml/10min
IS-IP828	LIAprim	Infusion pump 275ml, flow rate 6/8/10/12ml/h, PCA Bolus 2ml/15min