

Uprim
URINE BED & NIGHT BAGS
LEG BAGS
BODY WORN BAGS
PRODUCT CODES



Uprim PEDIATRIC URINE BAGS



Product code	Designation	Description
--------------	-------------	-------------

UR-BA01	Uprim	Paediatric, 100 mL
UR-BA02	Uprim	Paediatric, 200 ml

DESCRIPTION:

Urinary pediatric collection bag with safety adhesive, in 100 ml and 200 ml volume, opened sided face, for universal use (for girls and boys).

MATERIALS

PE, adhesive dressing

METHOD OF STERILISATION

Ethylene oxide

SHELF LIFE

5 years

Uprim LEG BAG / BODY WORN BAG



Product code	Designation	Description
UR-BA03	Uprim	Urine bag (three-chamber leg bag 500 ml, Flip type outlet)

DESCRIPTION:

500 ml, three-chamber, leg urine bag with Flip type outlet valve, tube O.D. 7 mm, length 45 cm. Made of white foil. Needleless adapter, connector with cap and anti-reflux valve.

MATERIALS

PVC, PE & PP

METHOD OF STERILISATION

Ethylene oxide

SHELF LIFE

5 years



Product code	Designation	Description
UR-BA04	Uprim	Urine bag (leg bag with Pull-Push outlet, 500 ml)
UR-BA05	Uprim	Urine bag (leg bag with Pull-Push outlet, 750 ml)

DESCRIPTION:

500 ml and 750 ml leg urine bag with Pull-Push outlet, tube O.D. 6.4 mm, length 30 cm. Made of white foil. Adapter with cap and anti-reflux valve.

MATERIALS

PVC, PE & PP

METHOD OF STERILISATION

Ethylene oxide

SHELF LIFE

5 years

Uprim LEG BAG / BODY WORN BAG



Product code	Designation	Description
UR-BA14	Uprim	Urine bag (three-chamber leg urine bag, 500 ml, T-valve outlet)

DESCRIPTION:

500 ml three-chamber leg urine bag with T-valve outlet, tube O.D. 7 mm, length 45 cm. Made of white foil. Adapter with cap and anti-reflux valve.

MATERIALS

PVC, PE & PP

METHOD OF STERILISATION

Ethylene oxide

SHELF LIFE

5 years

Uprim NIGHT AND DAY BAG



Product code	Designation	Description
UR-BA06	Uprim	Urine bag (night & day bag, 2000 ml Pull-Push outlet)
UR-BA07	Uprim	Urine bag (night & day bag, 2000 ml Screw type outlet)
UR-BA08	Uprim	Urine bag (night & day bag, 2000 ml T-valve outlet)

DESCRIPTION:

Large capacity bags (2000 ml) are recommended for post-operative patients (for patients confined to bed). This model has a urimeter for an accurate measurement of urine in intensive care patients. Large capacity bags are provided in white foil, with a O.D. 6.4 mm tube, length of 90 cm, cap adapter, anti-reflux valve and with different outlet taps.

MATERIALS

PVC, PE & PP

METHOD OF STERILISATION

Ethylene oxide

SHELF LIFE

5 years

Uprim NIGHT AND DAY BAG



Product code	Designation	Description
UR-BA10	Uprim	Urine bag (night & day bag, 1000 ml Pull-Push outlet)
UR-BA11	Uprim	Urine bag (night & day bag, 1000 ml Screw type outlet)
UR-BA12	Uprim	Urine bag (night & day bag, 1000 ml T-valve outlet)
UR-BA13	Uprim	Urine bag (night & day bag, 1500 ml without outlet)

DESCRIPTION:

Large capacity bags (1000 ml or 1500 ml) are recommended for post-operative patients (for patients confined to bed). The model has a urimeter for an accurate measurement of urine in intensive care patients. Large capacity bags are provided in white foil, with a O.D. 6.4 mm tube, length of 90 cm, cap adapter, anti-reflux valve and with different outlet taps.

MATERIALS

PVC, PE & PP

METHOD OF STERILISATION

Ethylene oxide

SHELF LIFE

5 years

Uprim URINARY DRAINAGE BAG



Product code	Designation	Description
UR-BA09	Uprim	Urine bag (drainage bag with handle, 2000 ml)

DESCRIPTION:

2000 ml urinary drainage bag with handle, air-vent and anti-reflux valve drip chamber, needleless sampling port with cap, 10 x 120 cm transparent inlet tube, slide clamp and bed sheet clamp, universal double hook and rope hanger, T-valve, rigid dog house.

MATERIALS

PVC, PE & PP, ABS

METHOD OF STERILISATION

Ethylene oxide

SHELF LIFE

5 years

Uprim URINE METER SYSTEMS



Product code	Designation	Description
--------------	-------------	-------------

UR-BA15	Uprim	Urine bag (urine meter system)
---------	-------	--------------------------------

DESCRIPTION:

500 ml + 1500 ml fourfold chamber, 500 cc urimeter with 1500 ml exchangeable urine bag, drip chamber with air-vent and anti-reflux valve. Needleless sampling port with cap. 10 * 150 cm transparent inlet tube. Slide clamp and bed sheet clamp. Universal double hook and rope hanger. T-valve outlet.

MATERIALS

PVC, PE & PP, ABS, Styrene-butadiene block copolymer, PA

METHOD OF STERILISATION

Ethylene oxide

SHELF LIFE

5 years



Product code	Designation	Description
--------------	-------------	-------------

UR-BA16	Uprim	Urine bag (urine meter system)
---------	-------	--------------------------------

DESCRIPTION:

400 ml + 2600 ml fourfold chamber, 400 cc urimeter with drip chamber 2600 ml urine bag. Air vent and anti-reflux valve. Needleless sampling port with cap. 10 * 150 cm transparent inlet tube. Slide clamp and bed sheet clamp. Universal double hook and rope hanger. T-valve outlet.

MATERIALS

PVC, PE & PP, ABS, Styrene-butadiene block copolymer, PA

METHOD OF STERILISATION

Ethylene oxide

SHELF LIFE

5 years