

VACUUM BLOOD COLLECTION TUBES PRODUCT CODES



VACUUM BLOOD COLLECTION TUBES

Tube type	Color code	Additive	Application
NO ADDITIVE TUBE		-	Immunology, serology, biochemical test
COAGULANT TUBE		Clot activator	
COAGULANT TUBE WITH GEL		Clot activator, Separating gel	
HEPARIN TUBE		Sodium Heparin	Biochemical tests
		Lithium Heparin	
		Lithium Heparin, Separating gel	
		Sodium Heparin, Separating gel	
PT TUBE		Sodium Citrate (1:9) 0.109mol/l or 0.129mol/l	Coagulation studies
GLUCOSE/LACTATE TUBE		Potassium oxalate, Sodium fluoride	Blood sugar examination
		EDTA, Sodium fluoride	
EDTA TUBE		EDTA K2	Blood routine examination, molecular diagnostics
		EDTA K3	
		EDTA K2, Separating gel	
ESR TUBE		Sodium Citrate (1:4) 0.109mol/l or 0.129mol/l	Blood cell sedimentation rate test

VACUUM BLOOD COLLECTION TUBES

Product code	Designation	Description
BC-TU01	VACUprim	Blood collection tube, 2ml, size 13x75mm, NO ADDITIVE
BC-TU02	VACUprim	Blood collection tube, 5ml, size 13x100mm, NO ADDITIVE
BC-TU21	VACUprim	Blood collection tube, 10ml, size 16x100mm, NO ADDITIVE
BC-TU03	VACUprim	Blood collection tube, 2ml, size 13x75mm, COAGULANT
BC-TU04	VACUprim	Blood collection tube, 5ml, size 13x100mm, COAGULANT
BC-TU22	VACUprim	Blood collection tube, 10ml, size 16x100mm, COAGULANT
BC-TU05	VACUprim	Blood collection tube, 2ml, size 13x75mm, GEL AND CLOT ACTIVATOR
BC-TU06	VACUprim	Blood collection tube, 5ml, size 13x100mm, GEL AND CLOT ACTIVATOR
BC-TU23	VACUprim	Blood collection tube, 10ml, size 16x100mm, GEL AND CLOT ACTIVATOR
BC-TU07	VACUprim	Blood collection tube, 2ml, size 13x75mm, EDTAK2
BC-TU08	VACUprim	Blood collection tube, 5ml, size 13x100mm, EDTAK2
BC-TU41	VACUprim	Blood collection tube, 5ml, size 13x100mm, EDTA K2 / SEPARATING GEL
BC-TU42	VACUprim	Blood collection tube, 10ml, size 16x100mm, EDTA K2 / SEPARATING GEL
BC-TU35	VACUprim	Blood collection tube, 2ml, size 13x75mm, EDTA K3
BC-TU36	VACUprim	Blood collection tube, 5ml, size 13x100mm, EDTA K3
BC-TU43	VACUprim	Blood collection tube, 5ml, size 13x100mm, EDTA K3 / SEPARATING GEL
BC-TU44	VACUprim	Blood collection tube, 10ml, size 16x100mm, EDTA K3 / SEPARATING GEL
BC-TU13	VACUprim	Blood collection tube, 2ml, size 13x75mm, LITHIUM HEPARIN
BC-TU14	VACUprim	Blood collection tube, 5ml, size 13x100mm, LITHIUM HEPARIN
BC-TU26	VACUprim	Blood collection tube, 10ml, size 16x100mm, LITHIUM HEPARIN
BC-TU16	VACUprim	Blood collection tube, 5ml, size 13x100mm, LITHIUM HEPARIN / SEPARATING GEL
BC-TU27	VACUprim	Blood collection tube, 10ml, size 16x100mm, LITHIUM HEPARIN / SEPARATING GEL
BC-TU17	VACUprim	Blood collection tube, 2ml, size 13x75mm, SODIUM HEPARIN
BC-TU18	VACUprim	Blood collection tube, 5ml, size 13x100mm, SODIUM HEPARIN
BC-TU28	VACUprim	Blood collection tube, 10ml, size 16x100mm, SODIUM HEPARIN
BC-TU39	VACUprim	Blood collection tube, 5ml, size 13x100mm, SODIUM HEPARIN/ SEPARATING GEL
BC-TU40	VACUprim	Blood collection tube, 10ml, size 16x100mm, SODIUM HEPARIN/ SEPARATING GEL

VACUUM BLOOD COLLECTION TUBES

Product code	Designation	Description
BC-TU19	VACUprim	Blood collection tube, 2ml, size 13x75mm, SODIUM CITRATE 1:4
BC-TU29	VACUprim	Blood collection tube, 3ml, size 13x75mm, SODIUM CITRATE 1:4
BC-TU45	VACUprim	Blood collection tube, 5ml, size 13x100mm, SODIUM CITRATE 1:4
BC-TU33	VACUprim	Blood collection tube, 2ml, size 13x75mm, double layer, SODIUM CITRATE 1:4
BC-TU34	VACUprim	Blood collection tube, 3ml, size 13x75mm, double layer, SODIUM CITRATE 1:4
BC-TU46	VACUprim	Blood collection tube, 5ml, size 13x100mm, double layer, SODIUM CITRATE 1:4
BC-TU20	VACUprim	Blood collection tube, 2ml, size 13x75mm, SODIUM CITRATE 1:9
BC-TU30	VACUprim	Blood collection tube, 3ml, size 13x75mm, SODIUM CITRATE 1:9
BC-TU31	VACUprim	Blood collection tube, 2ml, size 13x75mm, double layer, SODIUM CITRATE 1:9
BC-TU32	VACUprim	Blood collection tube, 3ml, size 13x75mm, double layer, SODIUM CITRATE 1:9