

CLEARspec VAGINAL SPECULA PRODUCT CODES



CLEARspec VAGINAL SPECULA

Product code	Designation	Description	Length (± 5mm)	Open position (± 5mm)	Width at tip (± 5mm)
GY-CS01	CLEARspec	Vaginal specula with middle screw (Type C) small.	100 mm	58 mm	23 mm
GY-CS02	CLEARspec	Vaginal specula with middle screw (Type C) medium.	110 mm	80 mm	30 mm
GY-CS03	CLEARspec	Vaginal specula with middle screw (Type C) medium long.	138 mm	97 mm	26 mm
GY-CS04	CLEARspec	Vaginal specula with middle screw (Type C) large.	117 mm	87 mm	35 mm
GY-CS05	CLEARspec	Vaginal specula with side screw (Type A) small.	108 mm	107 mm	20 mm
GY-CS06	CLEARspec	Vaginal specula with side screw (Type A) medium.	112 mm	108 mm	28 mm
GY-CS07	CLEARspec	Vaginal specula with side screw (Type A) large.	120 mm	109 mm	33 mm
GY-CS08	CLEARspec	Vaginal speculum tache type (Type B) small.	111 mm	79 mm	16 mm
GY-CS09	CLEARspec	Vaginal speculum tache type (Type B) medium.	113 mm	83 mm	27 mm
GY-CS10	CLEARspec	Vaginal speculum tache type (Type B) large.	120 mm	88 mm	32 mm
GY-CS11	CLEARspec	Vaginal speculum central pin (Type D) extra small.	95 mm	63 mm	15 mm
GY-CS12	CLEARspec	Vaginal speculum central pin (Type D) small.	106 mm	73 mm	20 mm
GY-CS13	CLEARspec	Vaginal speculum central pin (Type D) medium.	107 mm	73 mm	24 mm
GY-CS14	CLEARspec	Vaginal speculum central pin (Type D) large.	113 mm	80 mm	30 mm
GY-CS15	CLEARspec	Vaginal speculum trigger type (Type E) small.	115 mm	80 mm	22 mm
GY-CS16	CLEARspec	Vaginal speculum trigger type (Type E) medium.	117 mm	85 mm	27 mm
GY-CS17	CLEARspec	Vaginal speculum trigger type (Type E) large.	120 mm	89 mm	30 mm